

ZOLL AED 3™

The AED that leads the way

IMAGE



MESSAGING

SMART. READY. TRUSTED.

The AED that leads the way

The ZOLL AED 3™ defibrillator guides rescuers through the process of performing high-quality CPR and can deliver a potentially lifesaving shock to the heart when needed. With a ZOLL AED 3 and quick action, tragedy could be averted. ZOLL's proprietary Real CPR Help® technology helps guide rescuers in performing high-quality CPR in accordance with the most recent guidelines for depth and rate of compression. The innovative, intuitive design and enhanced features of the ZOLL AED 3 will give unexpected heroes both the confidence and knowledge needed to treat sudden cardiac arrest (SCA).

The ZOLL AED 3 defibrillator provides:

- Enhanced Real CPR Help, featuring a full-colour display with vivid rescue images, CPR cycle timer, and a large colour bar gauge that shows CPR compression depth.
- CPR Uni-padz™ electrodes that can be used on both adult and paediatric patients, so there is no need for two separate sets of pads.
- Cloud connectivity that can automatically report AED readiness.
- 5-year shelf life for pads and batteries increases readiness and decreases frequency of maintenance.

ZOLL AED 3™ Semi-automatic

Semi-automatic AED Designed for lay and community responders.

PN: 8501-001201-13

ZOLL AED 3™ Automatic

Fully automatic AED Designed for lay and community responders.

PN: 8501-001202-13


Inclusions

- Large colour touch screen display with vivid rescue images
- Enhanced Real CPR Help, CPR cycle timer, and a large colour bar gauge that shows CPR compression depth.
- RapidShock Analysis™ enables the industry's shortest rhythm analysis minimising pauses to life-saving CPR.
- CPR Uni-padz™ Universal (Adult/Paediatric) electrodes (PN#8900-000280 updated design), 5 Year Shelf Life.
- AED 3 Battery Pack ((PN#8000-000696)
- WiFi enabled for PlusTrac Cloud Connection (device readiness) and Case Review™ (clinical event data case push)
- One free year of PlusTrac, Operators Manual, Quick Start Guide and Rescue Reference Poster.
- 8 year standard warranty.

ZOLL AED 3™

The AED that leads the way

CONSUMABLES



[Dropbox link to Image](#)

CPR Uni-padz™ Universal electrodes (Adult/Paediatric)

PN: 8900-000280

- Universal Adult or Child Electrode Pads
- Chest sensor includes Real CPR Help real-time feedback for depth and rate of compressions during CPR
- 5-year shelf life
- Rescue Kit: includes gloves, scissors, razor, barrier mask, wet/dry wipe




[Dropbox link to Image](#)

ZOLL AED 3™ Battery Pack (Lithium Manganese Dioxide Battery Pack)

PN: 8000-000696

- Lithium battery pack for use with ZOLL AED 3
- Lasts up to five years
- Smart technology - Displays power level on LCD and can also report it to PlusTrac Program Management via Wi-Fi



ZOLL AED 3™

The AED that leads the way

OPTIONAL EXTRAS



[Dropbox link to Image](#)

AED 3™ Wall Mount Bracket (Device Only)

PN: 8000-001255

Dimensions: 45.7 cm W x 25.4 cm H x 35.6 cm D

Weight: 0.18 kg

- The Universal Wall Mount holds the ZOLL AED 3 alone (no bag)
- It is labeled with the green ILCOR AED symbol and includes mounting hardware, consumable tracking card, wall-mount template, and installation instructions



[Dropbox link to Image](#)

AED 3™ Case Wall Mount Bracket (Device Stored in Carry Case Only)

PN: 8000-001266

Dimensions: 45.7 cm W x 25.4 cm H x 35.6 cm D

Weight: 6.35 kg

- The Universal Wall Mount holds the ZOLL AED 3 in its carry case
- It is labeled with the green ILCOR AED symbol and includes mounting hardware, consumable tracking card, wall-mount template, and installation instructions.



[Dropbox link to Image](#)

Standard Surface Wall Cabinet

PN: 8000-001256

Dimensions: 48.9 cm W x 55.8 cm H x 22.8 cm D

Weight: 6.8 kg

- The Standard Surface Wall Cabinet is a battery operated, alarmed wall cabinet.
- Holds a ZOLL AED 3 plus a backup set of electrodes and a backup battery should you need them.
- Choose the surface-mounted model for easy installation.
- Includes an ILCOR (International Liaison Committee on Resuscitation) AED label set (French and English) and installation instructions.
- Wiring block in the strobe that allows connection of the cabinet to security system via low-voltage wiring.
- Strobe light sold separately (PN: 8000-001259)

ZOLL AED 3™

The AED that leads the way

OPTIONAL EXTRAS



[Dropbox link to Image](#)

ZOLL AED 3™ Carry Case

PN: 8000-001250

Dimensions: 33 cm W x 17.8 cm H x 30.5 cm D

Weight: 0.4 kg

- Holds ZOLL AED 3 or AED 3 BLS with pre-connected CPR Uni-padz plus a backup set of CPR Uni-pad.
- Keeps the AED 3 secure and accessible.
- Provides incremental protection from dust and water. Increases the AED's IP rating to IP67 when enclosed in the case.



[Dropbox link to Image](#)

NEW* ZOLL AED 3™ Soft Carry Case

PN: 8000-001280

- Supports quick access to the ZOLL AED 3 device without removing the case
 - No buttons, snaps or zippers. Simply lift the front cover to gain access to the AED and the USB port
- Additional insert on lower front flap for owner ID information or accessory expiration dates
- Adjustable shoulder strap
- Easy viewing of the status indicator
- Made of rugged material
 - Protects the ZOLL AED 3 from the elements
 - Easy to clean with soap/water, 1:10 solution bleach (one part bleach and nine parts water)
- Storage for accessories. Store as many as two sets of CPR Uni-Padz® (one pre-connected and one spare) or connect pads and spare battery

ZOLL AED 3™

The AED that leads the way

OPTIONAL EXTRAS



NEW* ZOLL AED Rescue Pack

PN: 8000-001468-01

Backpack Features:

- Main compartment features three removable stick-in-place pouches
- Internal Velcro lining strips to secure modules or other Velcro compatible items
- Side pockets for easy access to vital sign tools, badges, and smaller equipment
- Durable and cleanable Tarpaulin fabric, with an ounce urethane coating, for easier cleaning and protection
- A stow-away "Rain fly" is included for bad weather

Small Rigid Plastic Case

PN: 8000-001253

Dimensions: 35.7 cm W x 20.3 cm H x 45.7 cm D

Weight: 1.36 kg



- Holds ZOLL AED 3 or AED 3 BLS with pre-connected Uni-padz along with one spare battery.
- For lightweight maximum protection, choose the slim Hard-shell Carry Case.
- Keeps the AED 3 more secure and accessible.
- Provides incremental protection from dust and water. Increases the AED's IP rating to IP67 when enclosed in the case.

Large Rigid Plastic Case

PN: 8000-001254

Dimensions: 40.6 cm W x 20.3 cm H x 50.8 cm D

Weight: 1.8 kg



- Holds ZOLL AED 3 or AED 3 BLS with pre-connected Uni-padz along with a backup battery pack and CPR Uni-padz.
- Protect and carry everything needed for a vigorous rescue with the large Hard-shell Carry Case.
- Keeps the AED 3 most secure and accessible.
- Provides incremental protection from dust and water. Increases the AED's IP rating to IP67 when enclosed in the case.

ZOLL AED 3™ Trainer

The AED that leads the way

OPTIONAL EXTRAS



[Dropbox link to Image](#)

SMART. READY. TRUSTED.

The AED that leads the way

The ZOLL AED 3 Trainer is designed in a way to closely mimic the design of the clinical ZOLL AED 3, allowing a realistic training experience.

The AED 3 Trainer have Real CPR Help incorporated, allowing the trainee to receive real-time compression feedback during a training.

Features:

- Real CPR Help
- Pads on Manikin Detection
- Active LCD Display
- "On-trainer Configuration"
- USB Software Updates

The ZOLL AED 3™ Trainer

PN: 8028-000001-07

Inclusions

- ZOLL AED 3 Trainer CPR Uni-padz Electrode Training Harness
- Trainer Electrode Gels for ZOLL AED 3 Trainer CPR Uni-padz, one pair
- Trainer Electrode Liner for ZOLL AED 3 Trainer CPR Uni-padz, one pair
- Batteries (4 D-Cell)
- ZOLL AED 3 Trainer Operator's Guide
- 12 Months Factory Warranty



ZOLL AED 3™ Trainer

The AED that leads the way

CONSUMABLES




[Dropbox link to Image](#)

ZOLL AED 3 Trainer
CPR Uni-padz Electrode Training Harness
PN: 8028-000010



[Dropbox link to Image](#)

ZOLL AED 3 Trainer
CPR Uni-padz Electrode Gels, Pack of 5
PN: 8028-000014



[Dropbox link to Image](#)

Replacement Liner for ZOLL AED Trainer
CPR Uni-padz Electrode, each
PN: 8028-000013

ZOLL AED 3™

The AED that leads the way

LITERATURE & PRODUCT MANUALS



ZOLL AED DUAL BROCHURE

[Dropbox link to Image](#)



ZOLL AED 3 BROCHURE

[Dropbox link to Image](#)



ZOLL AED 3 ACCESSORIES AND TRAINING PRODUCTS

[Dropbox link to Image](#)



ZOLL AED 3 SPEC SHEET

[Dropbox link to Image](#)



ZOLL AED WIFI SET UP GUIDE

[Dropbox link to Image](#)



ZOLL AED 3 TRAINER SPEC SHEET

[Dropbox link to Image](#)

ZOLL AED 3™

The AED that leads the way

LITERATURE & PRODUCT MANUALS



ZOLL AED 3 QUICK START GUIDE

[Dropbox link to Image](#)



ZOLL AED 3 OPERATOR'S MANUAL

[Dropbox link to Image](#)



ZOLL AED 3 ADMINISTRATOR'S GUIDE

[Dropbox link to Image](#)



ZOLL AED 3 RESCUE REFERENCE GUIDE

[Dropbox link to Image](#)



ZOLL AED 3 AUTOMATIC RESCUE REFERENCE GUIDE

[Dropbox link to Image](#)



ZOLL AED 3 CPR POSTER ERC

[Dropbox link to Image](#)

ZOLL AED 3™

The AED that leads the way

VIDEO & WEBSITE LINKS



BUSINESS PARTNER RESOURCE PORTAL

Our new hub access product literature, sales tools, images and videos to help you sell our AEDs. We are continuously adding new tools and resources to help you succeed, so be sure to look out for your login details with this email and visit the portal frequently.

PORTAL access: portal.zoll.com/ps-apac/

Username: Use your work email address or sign up

Password: realcprhelp



ZOLL AED 3™ Promotional Video

Promotional video of the ZOLL AED 3 defibrillator showcasing ZOLL's;

- Leading edge AED technology
- Enhanced Real CPR Help System
- Key features and benefits (Paediatric capabilities, WiFi connectivity, advanced technology)


Please visit: info.zoll.com/aed-3



ZOLL AED 3™ Basic Set Up Videos

- ZOLL AED 3
- ZOLL AED 3 Basic Set-up
- ZOLL AED 3 Basic LCD Touchscreen Configuration

Please visit: bit.ly/2QaqqdY



Watch the ZOLL AED 3™ Training Videos

- ZOLL AED 3
- ZOLL AED 3 - Inside The Rescue (ERC)
- Real CPR Help Demo
- PlusTrac Video
- ZOLL AED 3 Basic Set-up
- ZOLL AED 3 Basic LCD Touchscreen Configuration

Please visit: bit.ly/2D7DVdb



ZOLLAED3.com

- ZOLL AED 3 Device Registration
- Documentation

Please visit: www.zollaed3.com/

For More Information To learn more about The ZOLL AED 3™, please contact your local representative or visit ZOLL online at bit.ly/2D5rpeg

SMART. READY. TRUSTED.
Designed for Unexpected Heroes™

ZOLL®

ZOLL®